



Fecal Dx* antigen testing

Detect parasite infections even before eggs are present.

IDEXX


Find what the microscope misses.¹⁻⁵

Forget about eggs. Fecal Dx* antigen testing detects antigens secreted directly from the infecting worm.

- + Identify infections prior to eggs being shed.
- + Avoid false-negative results caused by intermittent shedding of eggs.
- + Detect worms even if eggs are not present in the faecal sample.

Help prevent spread of disease to family members and pets.

Detecting infections sooner helps you treat effectively, reducing the risk of infection to other pets and family members.⁶



A single test finds hookworms, roundworms, whipworms and flea tapeworms.



Earlier detection than faecal flotation.¹⁻⁵




Detects 5X more infections than faecal flotation alone.^{2,3,5}


Set your sights on parasites.

Fecal Dx* antigen testing helps you detect more parasites earlier than ever before.

Detect more. Results indicate that adding Fecal Dx antigen testing to your faecal flotation testing may double the number of parasite infections detected.⁷



Detect earlier. Because of the lack of egg detection with faecal flotation testing during the prepatent period and single-sex infections, many parasite infections may go undetected for a period of time.¹



¹Data represents faecal testing performed during wellness visits.

Parasites are everywhere.

Pan-European DOGWALKS study confirms the need to test and treat even well-cared for dogs.⁸


- + Approximately 2470 fresh canine fecal specimens were collected.
- + 164 parks in 33 cities across 12 European countries were included.
- + Fecal Dx* antigen testing, *Giardia* antigen testing and faecal flotation by centrifugation were conducted.

Key results:

57%



57% of parks tested positive for an intestinal parasite.



Almost 1 in 4 dogs tested positive for intestinal parasites.

37.2%



28.7%



17.7%



Between 18% and 37% of parks tested positive for hookworms (37.2%), roundworms (28.7%) or whipworms (17.7%).



Test to protect with Fecal Dx antigen testing:



Fecal Dx* antigen testing catches infections missed by the microscope.¹⁻⁴




Parasite infections may go undetected, even in dogs receiving preventives.⁸



Ideally, dogs are protected by broad-spectrum, year-round prevention.

Only IDEXX Reference Laboratories offers a complete faecal testing solution:




Faecal flotation by centrifugation



Egg and parasite identification by expert technicians.



Fecal Dx antigen testing.



Common test codes for faecal testing.

For healthy dogs and cats:

**Fecal Dx* Antigen Panel (test code NEMAE;
addon: NEMAES)**

Hookworm, roundworm, whipworm and flea tapeworm
antigen immunoassays

For sick dogs and cats:

**Fecal Dx* Antigen Profile plus *Giardia*
(test code FECG)**

Giardia, hookworm, roundworm, whipworm and
leishmaniasis antigen immunoassays

**Fecal Dx* Antigen Panel plus, with *Giardia* and
leishmaniasis (test code FECGP)**

Faecal leishmaniasis, *Giardia*, hookworm, roundworm,
whipworm and leishmaniasis antigen immunoassays




Just start with VetConnect PLUS, your go-to diagnostic resource for efficiency and flexibility.

VetConnect* PLUS offers the diagnostic tools you need for fast, effective clinical decision-making.



Spot trends and abnormalities with at-a-glance patient historical results.

Easily compare results side-by-side or graph patient historical results for visual trend analyses.




Improve client communications and increase compliance.

Share diagnostic results and personalized client-friendly result summaries with one click.



Access differential diagnoses to help interpret results and determine next steps

Click to view guidance from IDEXX board-certified specialists using the latest industry protocols, research and publications.



Get results in one to three days for Fecal Dx* antigen testing.

Results are automatically added to the patient history in VetConnect PLUS.

To activate your VetConnect PLUS account, please complete the online registration form on idexx.se/vetconnectplus | idexx.dk/vetconnectplus | idexx.no/vetconnectplus | idexx.fi/vetconnectplus
Sign in at vetconnectplus.se | vetconnectplus.dk | vetconnectplus.no | vetconnectplus.fi

Manage patient diagnostics from your smartphone, tablet, browser or practice management system. Mobile app is available for iPhone* and Android* devices.



References

1. Elsemore DA, Geng J, Flynn L, Cruthers L, Lucio-Forster A, Bowman DD. Enzyme-linked immunosorbent assay for coproantigen detection of *Trichuris vulpis* in dogs. *J Vet Diagn Invest.* 2014;26(3):404–411.
2. Data on file at IDEXX Laboratories, Inc. Westbrook, Maine USA.
3. Data on file at IDEXX Laboratories, Inc. Westbrook, Maine USA. Aggregate detection of hookworm, roundworm, and whipworm infections.
4. Elsemore DA, Geng J, Cote J, Hanna R, Lucio-Forster A, Bowman DD. Enzyme-linked immunosorbent assays for coproantigen detection of *Ancylostoma caninum* and *Toxocara canis* in dogs and *Toxocara cati* in cats. *J Vet Diagn Invest.* 2017;29(5):645–653. doi:10.1177/1040638717706098
5. Elsemore D, Beall M, Bezold T, et al. Detection of *Dipylidium caninum* coproantigen in experimental and natural infections [AAVP Abstract 23]. Paper presented at: American Association of Veterinary Parasitologists 67th Annual Meeting; June 26, 2022; Snowbird, UT
6. CAPC guidelines: controlling internal and external parasites in U.S. dogs and cats. Companion Animal Parasite Council. Accessed March 9, 2023. www.petsandparasites.org/resources/capc-guidelines
7. Sweet S, Hegarty E, McCrann DJ, Coyne M, Kincaid D, Szlosek D. A 3-year retrospective analysis of canine intestinal parasites: fecal testing positivity by age, U.S. geographical region and reason for veterinary visit. *Parasit Vectors.* 2021;14(1):173. doi:10.1186/s13071-021-04678-6
8. Drake J, Sweet S, Baxendale K, et al. Detection of *Giardia* and helminths in Western Europe at local K9 (canine) sites (DOGWALKS Study). *Parasit Vectors.* 2022;15(1):311. doi:10.1186/s13071-022-05440-2

© 2023 IDEXX Laboratories, Inc. All rights reserved. • 09-2690380-01 • The IDEXX Privacy Policy is available at idexx.com.

*Fecal Dx, PetChek and VetConnect are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries. All other product and company names and logos are trademarks or registered trademarks of their respective holders.

IDEXX